



## User's Manual

# Human Meteorin Like ELISA Kit



DEIA-XY2224



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This Human Meteorin Like Protein (METRNL) ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural human Meteorin Like Protein from serum and plasma in a sandwich ELISA format.

### Storage

2-8°C

### Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

### Detection Range

0.977-125 ng/mL

### Sensitivity

150 pg/mL

### Specificity

Human Meteorin Like Protein

