



## User's Manual

# Mouse/Rat MFGE8 ELISA Kit

REF

DEIA-XY2225



96T

RUO



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This Rat/Mouse Milk fat globule-EGF factor 8 (MFGE8) ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural rat MFGE8 from milk, serum and plasma in a sandwich ELISA format.

### Storage

2-8°C

### Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

### Detection Range

312.5-20000 pg/mL

### Sensitivity

50 pg/mL

### Specificity

Rat, Mouse