



User's Manual

Human Total MMP1 ELISA Kit

REF

DEIA-XY2226



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of human total mmp-1 concentrations in cell culture supernates, serum and EDTA plasma.

Storage

2-8°C

Precision

Intra-assay Precision: 4-8%; Inter-assay Precision: 8-12%

Detection Range

156-20,000 pg/mL

Sensitivity

15 pg/mL

Specificity

Human Total MMP-1