



User's Manual

Human Total MMP9 ELISA Kit

REF DEIA-XY2230

 96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** info@creative-diagnostics.com  **Web:** www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative determination of human total mmp-9 concentrations in serum and EDTA plasma.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-12%

Detection Range

31.25-2000 pg/mL

Sensitivity

15.6 pg/mL

Specificity

Human MMP-9 (Pro- and active form)