



**User's Manual**

# Human Soluble Neuropilin-2 ELISA Kit



**DEIA-XY2231**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of human soluble neuropilin 2 concentrations in plasma and serum.

### Storage

2-8°C

### Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

### Detection Range

125-8000 pg/mL

### Sensitivity

50 pg/mL

### Specificity

Human soluble Neuropilin 2 only