



**User's Manual**

# Human Soluble Nidogen-2 ELISA Kit



DEIA-XY2232



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of human nidogen-2 concentrations in serum and plasma.

### Storage

2-8°C

### Precision

Intra-assay Precision: 6-8%; Inter-assay Precision: 8-12%

### Detection Range

3.12-200 ng/ml

### Sensitivity

50 pg/mL

### Specificity

Human Nidogen-2