



User's Manual

Human Osteocrin-prepro ELISA Kit

REF

DEIA-XY2233



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This Human Osteocrin-prepro ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural human Osteocrin-prepro from EDTA plasma in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

Detection Range

39-2500 ng/mL

Sensitivity

5 ng/mL