



User's Manual

Human Soluble RANK ELISA Kit

REF

DEIA-XY2237



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of human sRANK concentrations in cell culture supernates

Storage

2-8°C

Precision

Intra-assay Precision: 6-8%; Inter-assay Precision: 8-12%

Detection Range

62-4000 pg/mL

Sensitivity

31 pg/mL

Specificity

Human RANK only