



User's Manual

Human Soluble HSP70 ELISA Kit

REF

DEIA-XY2241



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of human shsp70 concentrations in cell lysates, serum and EDTA plasma.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

Detection Range

156.25-20000 pg/MI

Sensitivity

50 pg/mL

Specificity

Human, Mouse, Rat sHSP70