



**User's Manual**

# Human Soluble LOX-1 ELISA Kit

REF

DEIA-XY2244



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

This Human Soluble Lectin-like oxidized LDL Receptor 1 (LOX-1) ELISA Kit contains the necessary components required for the quantitative measurement of recombinant human soluble LOX-1 (human cells derived) and/or natural human soluble LOX-1 from serum and plasma samples in a sandwich ELISA format.

### Storage

2-8°C

### Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

### Detection Range

39-2500 pg/mL

### Sensitivity

7.8 pg/mL

### Specificity

Human Soluble LOX-1