



User's Manual

Methyltestosterone ELISA kit



DEIA-XY31



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA kit is used for the quantitative analysis of methyltestosterone in cattle urine or tissue samples.

Storage

2-8°C

Detection Limit

2.0 ppb for urine samples and 0.0625 ppb for tissue

Specificity

The kit uses a polyclonal rabbit antibody which shows the following cross-reactivities: methyltestosterone (100%), methandrolone (47%), β -testosterone (33%), stanozolol (18%), 19-nortestosterone (2.7%), dihydrotestosterone (2.3%), and β -trenbolone (<0.1%).



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com