



User's Manual

Sulfonamides II ELISA kit



DEIA-XY33



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The multi-screening sulfonamides II ELISA is a competitive immunoassay for tissue, milk, honey and urine.

Storage

2-8°C

Detection Limit

<15 ng/mL for the above mentioned matrices

Specificity

The antibody used in this kit has large cross-reactivities with 13 sulfonamides, i.e. sulfadiazine, sulfadoxine, sulfamethoxazole, sulfisoxazole, sulfapyridine, sulfachlorpyridazine, sulfathiazole, sulfamethizole, sulfamethoxypyridazine, sulfadimethoxine, sulfabenzamide, sulfamethoxydiazine and sulfamerazine.