



User's Manual

Nortestosterone ELISA kit



DEIA-XY35



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The 19-Nortestosterone ELISA is a competitive enzyme immunoassay based on antibodies directed against 17 β -19-Nortestosterone for milk and urine.

Storage

2-8°C

Detection Limit

0.5 ng/mL for urine or milk samples

Specificity

The ELISA kit contains an antibody which cross-reacts with 17 β -19-nortestosterone (100%), 17 α -19-nortestosterone (50%), nortestosterone-sulphate (100%), 19-nor-4-androstene-3,17-dione (240%), norgestrel (16%) and norethindrone (38%).