



## User's Manual

# Semicarbazide ELISA kit



DEIA-XY39



96T






This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The SEM ELISA is used for the quantitative analysis of SEM residues in tissue (muscle, liver), shrimps, egg (powder), milk, and urine.

### Storage

2-8°C

### Detection Limit

<0.3 ng/g for the above mentioned matrices

### Specificity

The kit uses a monoclonal mouse antibody which shows the following cross-reactivities: SEM (100%), AMOZ (<0.01%), AOZ (<0.01%), and AHD (<0.01%).

