



User's Manual

Stanozolol ELISA kit



DEIA-XY40



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The stanozolol ELISA is used for the quantitative analysis of stanozolol in cattle urine or faeces samples.

Storage

2-8°C

Detection Limit

1 CC β for urine or faece samples

Specificity

The kit uses a polyclonal rabbit antibody which shows the following cross-reactivities: stanozolol (100%), 16 β -hydroxystanozolol (100%), 4 β -hydroxystanozolol (0.3%), 3'-hydroxystanozolol (<0.01%), methyltestosterone (<0.01%), and testosterone (<0.01%).