



User's Manual

Human APL1 beta 27 ELISA kit



DEIA-XY67



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This assay kit is capable for the quantitative determination human APL1 β 27 in cerebrospinal fluid and cell culture supernatant.

Storage

2-8°C

Precision

Intra-Assay Precision: 6.6%-7.1%

Inter-Assay Precision: 4.8%-8.4%

Sensitivity

2.72 pg/mL

Specificity

Human APL1 β 27: 100 %

Human APL1 β 25: <0.1

Human APL1 β 28: <0.2

