



User's Manual

Bordetella FHA IgA ELISA kit



DEIA-XY69



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Enzyme immunoassay for the quantitative determination of IgA antibodies against Bordetella pertussis FHA (filamentous hemagglutinin) in human serum and plasma.

Storage

The kit is shipped at ambient temperature and should be stored at 2-8°C. Keep away from heat or direct sun light.

Precision

Intra-Assay Precision: 6.0%-12.6%

Inter-Assay Precision: 8.1%-17.1%

Sensitivity

100%

Specificity

100%