



User's Manual

Rabbit CRP ELISA kit

REF

DEIA-XY76



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Enzyme Immunoassay for the quantitative determination of rabbit c-reactive protein (CRP) in serum.

Storage

1. The lyophilized reference standard should be stored at minus 20°C for optimum stability (it can be safely shipped at 2-8°C).
2. The remainder of the kit should be stored at 2-8°C and the microtiter plate should be kept in a sealed bag with desiccant to minimize exposure to damp air.
3. Test kits will remain stable for six months from the date of purchase, provided that the components are stored as described above.