



User's Manual

Phenylalanine neonate (with Membrane Plate) ELISA kit



DEIA-XY8



192T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Phenylalanine Microwell Enzyme Assay is designed for the quantitative determination of phenylalanine in neonatal blood spots.

It is for use as an aid in screening newborns for elevated phenylalanine levels.

Storage

2-8°C