



User's Manual

ADAM12 ELISA Kit



DEIA-XY94



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Enzyme-linked immunosorbent assay (ELISA) for the determination of ADAM12 in cell culture supernatants, serum, plasma, or urine.

Storage

2-8°C

Detection Range

156-10000 pg/mL

Sensitivity

The minimum detectable dose of this ADAM12 ELISA typically is less than 10 pg/mL.

Specificity

This ADAM12 ELISA allows for the precise measurement of endogenous levels of natural and/or recombinant ADAM12 proteins within the range of 156 pg/mL-10000 pg/mL.