



User's Manual

Human ANCA Screen ELISA Kit



DEIA-XY98



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The ANCA screen assay is an enzyme immunoassay (EIA) intended to screen for the presence of IgG class autoantibodies against PR3 and MPO in human serum or plasma.

Storage

2-8°C

Specificity

The microplate is coated with a mixture of PR3 and MPO antigens, highly purified by affinity chromatography. The ANCA screen test is specific only for autoantibodies directed to these antigens. No cross-reactivities have been observed.