



User's Manual

CALR ELISA Kit



DEIA-XYA1003



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The CALR Cell-Based ELISA Kit is a convenient, lysate-free, high throughput and sensitive assay kit that can monitor CALR protein expression profile in cells. The kit can be used for measuring the relative amounts of CALR in cultured cells as well as screening for the effects that various treatments, inhibitors (ie. siRNA or chemicals), or activators have on CALR.

Storage

4°C/6 Months