



User's Manual

JunB ELISA Kit



DEIA-XYA1083



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION


Intended Use

The JunB DNA-Binding ELISA Kit is intended for qualitative analysis of activated transcription factor JunB in nuclear and cell lysates.

Storage

4°C/6 Months

Reference Values



Specificity

The JunB DNA-Binding ELISA Kit detects endogenous levels of total JunB protein.