



User's Manual

Keratin 16 ELISA Kit



DEIA-XYA1096



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Keratin 16 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of Keratin 16 and the effects that certain stimulation conditions have on Keratin 16 expression in different cell lines.

Storage

4°C/6 Months