



User's Manual

Keratin 18 ELISA Kit



DEIA-XYA1098



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Keratin 18 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of Keratin 18 and the effects that certain stimulation conditions have on Keratin 18 expression in different cell lines.

Storage

4°C/6 Months