



User's Manual

KLF7 ELISA Kit



DEIA-XYA1114



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION


Intended Use

The KLF7 DNA-Binding ELISA Kit contains components necessary for qualitative detection of active transcription factor KLF7 in eukaryotic nuclear or cell lysates.

Storage

4°C/6 Months

Reference Values



Specificity

The KLF7 DNA-Binding ELISA detects endogenous levels of KLF7.