



User's Manual

KSR ELISA Kit

REF DEIA-XYA1117

Σ 96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** info@creative-diagnostics.com  **Web:** www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The KSR Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of KSR and the effects that certain stimulation conditions have on KSR expression in different cell lines.

Storage

4°C/6 Months