



User's Manual

LKB1 ELISA Kit



DEIA-XYA1138



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** info@creative-diagnostics.com  **Web:** www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The LKB1 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of LKB1 and the effects that certain stimulation conditions have on LKB1 expression in different cell lines.

Storage

4°C/6 Months