



User's Manual

## MAP3K7 ELISA Kit

**REF** DEIA-XYA1160

 96T


**RUO**



This product is for research use only and is not intended for diagnostic use.



For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The MAP3K7 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of MAP3K7 and the effects that certain stimulation conditions have on MAP3K7 expression in different cell lines.

### Storage

4°C/6 Months