



## User's Manual

# MEK1/2 ELISA Kit



**DEIA-XYA1198**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The MEK1/2 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of MEK1/2 and the effects that certain stimulation conditions have on MEK1/2 expression in different cell lines.

### Storage

4°C/6 Months