



User's Manual

MTA1 ELISA Kit

REF DEIA-XYA1238


 96T



RUO



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The MTA1 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of MTA1 and the effects that certain stimulation conditions have on MTA1 expression in different cell lines.

Storage

4°C/6 Months