



User's Manual

Mouse Eotaxin (Luminometer) ELISA Kit

REF

DEIA-XYA1244



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Eotaxin (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse Eotaxin concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Specificity

The Mouse Eotaxin ELISA is capable of recognizing both recombinant and naturally produced Mouse Eotaxin proteins. The antigens listed below were tested at 50 ng/mL and exhibited less than 1% cross reactivity.

Human: Eotaxin The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: Eotaxin-2, Eotaxin-3, MCP-1, MCP-2, MCP-3, MCP-4

Mouse: C10, CTACK, Eotaxin-2, Exodus-2, JE, MCP-2, MCP-5, MDC, MEC, MIP-1alpha, MIP-1beta, MIP-1gamma, MIP-3beta, RANTES

Rat: MCP-1