



User's Manual

Mouse G-CSF (Colorimetric) ELISA Kit

REF

DEIA-XYA1247



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The G-CSF (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse G-CSF concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Detection Range

32-2000 pg/mL

Specificity

The Mouse G-CSF ELISA is capable of recognizing both recombinant and naturally produced Mouse G-CSF proteins. The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: CT-1, G-CSF, GM-CSF, IL-3, IL-6, Leptin, M-CSF, SCF

Mouse: CT-1, GM-CSF, IL-3, IL-6, Leptin, M-CSF, SCF

Rat: GM-CSF, IL-3beta, IL-6, Leptin