



**User's Manual**

# Mouse MCP-3 (Colorimetric) ELISA Kit



**DEIA-XYA1261**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The MCP-3 (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse MCP-3 concentrations within any experimental sample including cell lysates, serum and plasma.

### Storage

4°C/6 Months

### Detection Range

8-2000 pg/mL

### Specificity

The Mouse MCP-3 ELISA is capable of recognizing both recombinant and naturally produced Mouse MCP-3 proteins. The antigens listed below were tested at 50 ng/mL and exhibited less than 1% cross reactivity.

Mouse: MDC The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: MCP-1, MCP-2, MCP-3, MCP-4

Mouse: C-10, CTACK, MDC-2, Exodus-2, JE, MCP-2, MCP-5, MDC, MIP-1alpha, MIP-1beta, MIP-1gamma, MIP-3alpha, MIP-3beta, RANTES

Rat: MCP-1

