



User's Manual

Mouse MIP-1beta (Luminometer) ELISA Kit



DEIA-XYA1266



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The MIP-1beta (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse MIP-1beta concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Specificity

The Mouse MIP-1beta ELISA is capable of recognizing both recombinant and naturally produced Mouse MIP-1beta proteins. The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: MIP-1alpha, MIP-1beta, MIP-3, MIP-3alpha, MIP-3beta, MIP-4, MIP-5

Mouse: MIP-1alpha, MIP-1gamma, MIP-2, MIP-3alpha, MIP-3beta

Rat: MIP-1alpha, MIP-1beta