



User's Manual

Mouse MIP-1beta (Colorimetric) ELISA Kit



DEIA-XYA1267



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The MIP-1beta (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse MIP-1beta concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Detection Range

4-500 pg/mL

Specificity

The Mouse MIP-1beta ELISA is capable of recognizing both recombinant and naturally produced Mouse MIP-1beta proteins. The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: MIP-1alpha, MIP-1beta, MIP-3, MIP-3alpha, MIP-3beta, MIP-4, MIP-5

Mouse: MIP-1alpha, MIP-1gamma, MIP-2, MIP-3alpha, MIP-3beta

Rat: MIP-1alpha, MIP-1beta