



**User's Manual**

# Mouse MIP-2 (Luminometer) ELISA Kit



**DEIA-XYA1268**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The MIP-2 (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse MIP-2 concentrations within any experimental sample including cell lysates, serum and plasma.

### Storage

4°C/6 Months

### Specificity

The Mouse MIP-2 ELISA is capable of recognizing both recombinant and naturally produced Mouse MIP-2 proteins. The antigens listed below were tested at 50 ng/mL and exhibited less than 1% cross reactivity.

Human: GRObeta

Mouse: I-TAC

Rat: GRObeta

The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: GRO, GROgamma, NAP-2

Mouse: CXCL-16, IP-10, JE, KC, LIX, MIG, MIP-1alpha, MIP-1beta, MIP-1gamma, MIP-3alpha, MIP-3beta, SDF-1alpha, SDF-1beta

Rat: GRO/KC