



User's Manual

Mouse MIP-2 (Colorimetric) ELISA Kit



DEIA-XYA1269



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**



Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com



Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The MIP-2 (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse MIP-2 concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Detection Range

32-2000 pg/mL

Specificity

The Mouse MIP-2 ELISA is capable of recognizing both recombinant and naturally produced Mouse MIP-2 proteins. The antigens listed below were tested at 50 ng/mL and exhibited less than 1% cross reactivity.

Human: GRObeta

Mouse: I-TAC

Rat: GRObeta The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: GRO, GROgamma, NAP-2

Mouse: CXCL-16, IP-10, JE, KC, LIX, MIG, MIP-1alpha, MIP-1beta, MIP-1gamma, MIP-3alpha, MIP-3beta, SDF-1alpha, SDF-1beta

Rat: GRO/KC