



User's Manual

Mouse RANTES (Colorimetric) ELISA Kit

REF

DEIA-XYA1271



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The RANTES (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse RANTES concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Detection Range

16-2000 pg/mL

Specificity

The Mouse RANTES ELISA is capable of recognizing both recombinant and naturally produced Mouse RANTES proteins. The antigens listed below were tested at 50 ng/mL and exhibited 100% cross reactivity.

Rat: RANTES

The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: MEC, RANTES, TECK

Mouse: C-10, RANTES, JE, MCP-2, MCP-3, MCP-5, MEC, MIP-1alpha, MIP-1beta, MIP-1gamma, MIP-2, MIP-3alpha

Rat: MCP-1, MIP-1alpha, MIP-1beta