



User's Manual

Mouse SCF (Luminometer) ELISA Kit



DEIA-XYA1272



96T




This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The SCF (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse SCF concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Specificity

The Mouse SCF ELISA is capable of recognizing both recombinant and naturally produced Mouse SCF proteins. The antigens listed below were tested at 50 ng/mL and exhibited 100% cross reactivity.

Mouse: SCF

Rat: SCF

The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: flt3-L, GM-CSF, IL-1alpha, IL-1beta, M-CSF

Mouse: IL-1alpha, IL-1beta, IL-4 Rat: IL-1alpha, IL-1beta