



User's Manual

Mouse sRANK Ligand (Luminometer) ELISA Kit

REF

DEIA-XYA1274



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The sRANK Ligand (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse sRANK Ligand concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Specificity

The Mouse sRANK-Ligand ELISA is capable of recognizing both recombinant and naturally produced Mouse sRANK-Ligand proteins. The antigens listed below were tested at 50 ng/mL and exhibited 100% cross reactivity.

Rat: sRANKL

The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: sFas, sFasL, sCD40L, sRANKL, sRANKR, TNFalpha, TNFbeta, sTRAIL

Mouse: April, sCD40L, LIGHT, sRANKL, TNFalpha