



User's Manual

Mouse VEGF (Colorimetric) ELISA Kit



DEIA-XYA1281



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The VEGF (Mouse) ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Mouse VEGF concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Detection Range

12-5000 pg/mL

Specificity

The Mouse VEGF ELISA is capable of recognizing both recombinant and naturally produced Mouse VEGF proteins. The antigens listed below were tested at 50 ng/mL and exhibited 100% cross reactivity.

Human: VEGF121, VEGF165

The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: EGF, EG-VEGF, flt3L, GM-CSF, PDGF-AA, PDGF-AB, PDGF-BB, PIGF, SCF

Mouse: EGF, flt3L, GM-CSF, SCF

Rat: SCF