



User's Manual

NCOR2 ELISA Kit



DEIA-XYA1306



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION


Intended Use

The MNCOR2 DNA-Binding ELISA Kit contains components necessary for qualitative detection of active transcription factor NCOR2 in eukaryotic nuclear or cell lysates.

Storage

4°C/6 Months

Reference Values



Specificity

The NCOR2 DNA-Binding ELISA Kit detects endogenous levels of total NCOR2 protein.