



User's Manual

NKX6.3 ELISA Kit



DEIA-XYA1328



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The NKX6.3 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of NKX6.3 and the effects that certain stimulation conditions have on NKX6.3 expression in different cell lines.

Storage

4°C/6 Months