



User's Manual

NR1I3 ELISA Kit



DEIA-XYA1343



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION


Intended Use

The NR1I3 DNA-Binding ELISA Kit contains components necessary for qualitative detection of active transcription factor NR1I3 in eukaryotic nuclear or cell lysates.

Storage

4°C/6 Months

Reference Values



Specificity

The NR1I3 DNA-Binding ELISA detects endogenous levels of total NR1I3 protein.