



User's Manual

Nuclear Receptor NR4A1 ELISA Kit



DEIA-XYA1350



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Nuclear Receptor NR4A1 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of Nuclear Receptor NR4A1 and the effects that certain stimulation conditions have on Nuclear Receptor NR4A1 expression in different cell lines.

Storage

4°C/6 Months