



## User's Manual

# OCT3 ELISA Kit



**DEIA-XYA1359**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION


### Intended Use

The OCT3 DNA-Binding ELISA Kit contains components necessary for qualitative detection of active transcription factor OCT3 in eukaryotic nuclear or cell lysates.

### Storage

4°C/6 Months

### Reference Values



### Specificity

The OCT3 DNA-Binding ELISA detects endogenous levels of total OCT3 protein.