



## User's Manual

# p63 ELISA Kit



DEIA-XYA1382



96T





This product is for research use only and is not intended for diagnostic use.


For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The p63 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of p63 and the effects that certain stimulation conditions have on p63 expression in different cell lines.

### Storage

4°C/6 Months