



## User's Manual

# p70 S6 Kinase ELISA Kit

REF

DEIA-XYA1386



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The p70 S6 Kinase Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of p70 S6 Kinase and the effects that certain stimulation conditions have on p70 S6 Kinase expression in different cell lines.

### Storage

4°C/6 Months