



User's Manual

PDE4D ELISA Kit



DEIA-XYA1422




96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The PDE4D Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of PDE4D and the effects that certain stimulation conditions have on PDE4D expression in different cell lines.

Storage

4°C/6 Months